

Hi Glow!

Thermal Cold Beauty
Massager



USER MANUAL

DESCRIPTION OF THE DEVICE


1. ION+: Using the absorption power of positive ions, it removes dirt and aging skin behind the pores. It has higher detergency than normal makeup remover and absorbs clogged pores and dark spots that cannot be removed with a hand-washed face, leading to freshly born skin.

2. ION-: By passing a weak current, vitamin C and water-soluble beauty ingredients are ionized and delivered to the back of the skin (dermis layer), leading to a more transparent and moist skin.

3. Sonic vibration massage: It can be expected to activate facial muscle cells to improve metabolism and reduce skin swelling.

4. LED light beauty treatment:
*Red Light: It penetrates deep into the skin, promotes microcirculation through photochemistry, and stimulates metabolism and rejuvenation,Also promotes the growth of collagen.
*Blue Light: Acne is effectively matched to the absorption peak of acne propionic acidaceae, penetrating 0.25mm from the surface of the skin, and the cell disinfection function begins to act, destroying and sterilizing the bacterial habitat.
*Green Light: Calm and repair, mainly used for soothing sensitive skin.

DEVICE STRUCTURE



①The Heating and Cooling Panel

②LED Beads

③LCD Display (Four Gears Caption)

④Power Botton (Gears Adjustment)

⑤Charging Jack

⑥Device Holder(Base)

MODE AND DISPLAY

1. Long press the power button for 3 seconds to turn on/off the device. **Four gears** can be selected. Press the power button to switch the gears.

2. After turn on the device , **“WELCOME”** is displayed;

3. Gear 1: Heating (temperature is 35-39 ℃)+ Vibration, **“DEEP CLEAN”** is displayed;

4. Gear 2: Heating(temperature is 40-45 ℃)+ Red light + Vibration, **“NUTRITION”** is displayed;

5. Gear 3: Cooling(temperature is 12-15 ℃)+ Blue light, **“BOOST”** is displayed;

6. Gear 4: Vibration + Green light, **“VIBRATION”**is displayed;

7. Low Battery: LED light flashing (Depends on the current mode) , Power off after 10 seconds.

Charging: Red light flashing on the back of the device

Full Charged: Solid blue light on the back of the device

NOTE: Gear 1&2 will intelligent automatically shut off after 8 minutes use.
Gear 3 will intelligent automatically shut off after 3 minutes use. if turn on device but not use it in 12 minutes ,the device will shut off too. All are programmed for protect the battery.

HOW TO USE

Step 1 : After cleansing, apply the gel,lotion or essence on face.



Step 2 : Press the power button 3 seconds to turn on the device, select the appropriate gear, and use dot matrix pressing method or other face massage methods on our skin.(following use instructions for reference)
Gear 1 & 2 for promoting the absorption of skin care products.
Gear 3 & 4 for shrinking pores,firming skin, and effective skin care.

USE INSTRUCTIONS

**Forehead**
Glide the device from ajina to temple gently and slowly

**Nose**
Glide the device from nose bridge to the mouth direction

**Cheek**
Glide ttom the device from the bottom of the cheek to the top side

**Chin**
Glide the device by following the chin outline

**Tarea**
Glide the device by following the directions showed in the picture

**Lymphatic detoxification**
Glide the device from the both sides of face, then follow the directions of neck in top-down direction

PRODUCT PARAMETES

Product Name	Thermal Cold Beauty Massager
Material	ABS + Aluminum Alloy
Vibration Frequency	8000RPM
Heating Temperature	35-45 ℃
Cooling Temperature	12-15 ℃
Battery Capacity	450mA
Product Weight	160g
Product Dimension	140*52*28mm
Voltage	3.7v-5v
Charging Time	2-3 hours
Warranty	1 year
Technology	Ion+- & Heating & Cooling & LED light therapy (Red+blue+green) &Ultrasonic Vibration
Wavelength	Red Light: 622nm ± 10nm Blue Light: 465nm ± 10nm Green Light: 530nm ± 10nm

ATTENTION

1. Please keep it out of the reach of infants.

2. Please do not use following symptoms: heart disease patients. Severe patients, acute patients, infectious patients, allergic physical patients, dermatitis,etc.

3. Please do not use during pregnancy and lactation.

4. Please do not use in the following areas: wound, plastic surgery, eyeball.

5. Keep away from fire and high temperature.

6. Please stop using the device when it is abnormal.

WARNING

1. Please do not use other charging cables to avoid damaging the device.

2. Please do not use disassemble or repair this product by yourself to avoid damage.

3. Please avoid falling and hitting the product to avoid damage.

4. Do not rinse the product directly into water.

CLEANING AND MAINTENANCE

1. Wipe it off after use immediately.

2. This product can not be washed with water directly.

3. Do not use the sharp or corrosive or volatile diluent to clean the product.

ENVIRONMENTAL PROTECTION INSTRUCTIONS

1. This product is in the accordance with the regulations and of the management method for the restricted use of hazardous substances in the electrical and electronic products.

2. Consumers will not have the leakage and precipitation of harmful substances during normal use, consumer can be used safely.

3. The battery should be disposed of safely.

DISCLAIMER

All users of this device do so at their own risk.The manufacturer reserves the right to revise this document and may make changes on any occasion in the contents there of without the requirement to notify any person of such revision or changes. Model may be altered without notice.