

✱ Sales Points ✱

Electric Vibration LED Light Therapy Hot & Cold Nutrition Lead-in Device Ion Face Massager



① **Korean creativity, hand-drawn product prototypes, 100% original.** New and unique Appearance **PATENTED** Product for 2022 market. Just come to market November 2021. (CE & ROHS Certified)

2 , **Breakthrough 1:** Private model creation, only for light and thin, ultra-thin full-effect body, original grip, simple and generous appearance, high recognition, and can attract people at first glance.

Breakthrough 2: Innovatively equipped with cooling, heating, vibration, and phototherapy functions in the thin and light body , it is the one and only in the market. we do not follow the mess trend, Not tacky, no styling Collision



③**Material:** environmental ABS(shell) +14pcs Lamp beads + Aluminum Alloy (head)-


Higher Grade, No Corrode , No Oxidize, No Rust.Safe to our skin.



④ **Design:** **Portable and Luxury, Korea designed , Minimalist gold and silver with a white body, simple without losing personality, not easy to scratch,(looks like a Telephone)**

With a base; helps device stand steady~very convenient.

5. One-key operation, with built-in LCD display screen, displaying gears from inside to outside, making the operation layout modern and waterproof. The four corners of the



fuselage fit together without any gaps or looseness.

The head is safely raised a distance to better fit the skin and deliver nutrients point-to-point.

6, Ultra-fast heating and cooling control and gear switching, 4-speed cold and heat lifting, 3D vibration massage introduction. Subvert your perception and experience.

- * 35-39 °C warm skin to help the skin remove dirt
- * 40-45 °C Maintain warmth to help introduce nourishment
- * 12-15 °C icy feeling cools down and tightens skin
- * 12000RPM vibration massage to relieve skin fatigue

7, Product stability: Circuit temperature sensor (thermistor) is used as an additional protective measure to prevent excessive head temperature and low-temperature combustion in warm mode

②Specifications:

Product Name: Hot & Cold Nutrition Lead-in Device

Product Weight: Appr. 160g

Product Weight: Appr. 160g

Product Dimension: Appr. 140*52*28mm

Material: ABS+Zinc alloy

Power: 0.25w-2.5w

Battery: 450mA

Charging Time: Appr. 3 Hours

Technology: Ion+- & Heating & Cooling & LED light therapy (Red+blue+green) &Ultrasonic Vibration

Voltage: 3.7v-5v

Adapter requirements: 5V/1A

Working current: 0.05a-1.2a

Charging: USB rechargeable (Type C)

③Features:

<1>Vibration: 8000RPM/Min

<2>  Heating: 35-45°C

<3>  Cooling: 12-15°C

<4> Led light therapy:



{ **Red light: wavelength 622nm ± 10nm** }



{ **Blue light: wavelength 465nm ± 10nm** }



{ **Green light: wavelength 530nm ± 10nm** }

④Modes+Function:

Level 1: WELCOME (Turn On and Start Work) Funtion: /

Level 2: DEEP CLEAN <1>Ion+ Export*Ion- Import

<2>Warm (35-39°C) +vibration

(Turn off after using 8mins automatically)

Function: Ultrasonic cleaning; Help relieve muscles and promote pore absorption.

Level 3: NUTRITION Heating(40-45°C)+Red LED light+Vibration(Turn off after using 8mins automatically)

Function: Deep absorption skin care products, Anti-aging & wrinkle reducing by stimulating fibroblast cells to produce & secrete collagen.

Level 4: BOOST Cooling + Blue LED light (Turn off after using 8mins automatically)

Function: Shrink Pores and keep skin moist, Prevents and treats mild to moderate acne by killing acne causing bacteria ;Also controls sebum/oil production.

NOTE: Compared with other devices with cooling function, this device will get very cold once turn on this mode. Cooling speed is very fast ! ❄️

Level 5: VIBRATION Vibration+Green LED light

Function: Helps reduce hyperpigmentation; Help relieve muscles and promote pore absorption.

Timing Function to protect product :

It will automatically shut down if it is not used after two minutes !

Functions with good combination. **Small figure, big energy.** this device will improve fine lines, lift, tone, tighten, and provide lymphatic drainage to give an overall radiant glowing skin look.



⑤Item description & Product Details:

✓ ION Introduction

By passing a weak current, vitamin C and water-soluble beauty ingredients are ionized and delivered to the back of the skin (dermis layer), leading to a more transparent and moist skin.

✓ ION Derivation

Using the absorption power of positive ions, it removes dirt and aging skin behind the pores. It has higher detergency than normal makeup remover and absorbs clogged pores and dark spots that cannot be removed with a hand-washed face, leading to freshly born skin.

✓ Sonic vibration massage

It can be expected to activate facial muscle cells to improve metabolism and reduce skin swelling.

✓ RED light beauty treatment (wavelength 622nm $\pm$ 10nm)

It penetrates deep into the skin, promotes microcirculation through photochemistry, and stimulates metabolism and rejuvenation. Red light irradiation effectively activates collagen, quickly removes waste and dead cells from cell tissues and stimulates the secretion of fibroblast growth factor, so red light improves skin damage and fine lines. , The elasticity returns to the skin, the skin becomes dull and soft, moist and smooth.

✓ BLUE light beauty treatment (wavelength 465nm $\pm$ 10nm)

Acne is effectively matched to the absorption peak of acne propionic acidaceae, penetrating 0.25mm from the surface of the skin, and the cell disinfection function begins to act, destroying and sterilizing the bacterial habitat. As well as removing from the skin, you can expect to eliminate scars.

✓ Green light beauty treatment (wavelength 5305nm $\pm$ 10nm)

Helps reduce hyperpigmentation; Help relieve muscles and promote pore absorption.

✓ HOT function (35-45 °C, temperature will change hotter)

Cleansing and moisturizing care can be done comfortably while warming the skin.

✓ COLD function (12-15 °C temperature will change colder)

It cools your skin with a cool cooling effect.

✓ **Four-mode authentic aging care at home!**

⑥ Package included :

1×Device

1×Type-C Cable

1×Instruction Manual

1×Box

⑦ Customize service:

Logo on device customize (MOQ 50-100PCS)

Package customize (MOQ 100PCS-500PCS)

Instruction customize (MOQ 200PCS)

Other accessories customize (MOQ negotiable)

Device color customize (MOQ 1000PCS)

⑧ Details:

